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Distributed By
Medela Inc. USA
1101 Corporate Drive
McHenry, IL 60050
USA
Phone: +1 877 735 1626 Fax+1 800 995 7867
info-healthcare@medela.com
www.medela-healthcare.us

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Invia® – Silverlon NPWT Antimicrobial Wound Contact Dressing

Device Description

Invia – Silverlon NPWT Antimicrobial Wound Contact Dressing is a single layer of knitted nylon fiber substrate coated with metallic silver. This Antimicrobial Wound Contact Dressing model is knitted to a specified knit pattern designed to allow optimal flow rate of wound exudates when used as the primary wound contact layer in a NPWT dressing set. It is composed of 100% multifilament 15 denier nylon fiber, which is circumferentially coated with metallic silver at the rate of 0.68 mg/cm² (6.8 g/m²).

Ref. No	Size
087.7001	4" x 5" (10 x 12 cm)
087.7002	5" x 8" (12 x 20 cm)

Invia – Silverlon NPWT Antimicrobial Wound Contact Dressings deliver antimicrobial silver ions in the dressing when activated by moisture. The silver ions in the dressing kill wound bacteria held in the dressing and provide an antimicrobial barrier for bacterial penetration of the dressing which may help reduce infection. Invia - Silverlon NPWT Antimicrobial Wound Contact Dressings have been tested in vitro and found effective against following microorganisms:

- Gram Positive Bacteria - Methicillin Resistant Staphylococcus aureus (MRSA), Vancomycin Resistant Enterococcus (VRE), Staphylococcus epidermidis
- Gram Negative Bacteria - Escherichia coli (EColi), Pseudomonas aeruginosa, Enterobacter cloacae
- Yeast / Mold - Candida albicans, and Aspergillus brasiliensis

Invia – Silverlon NPWT Antimicrobial Wound Contact Dressings may be used with all Invia types of negative pressure wound therapy foam/gauze dressings and equipment.

Invia – Silverlon NPWT Antimicrobial Wound Contact Dressings have been subjected to independent standard in vitro and in vivo biocompatibility tests; cytotoxicity, sensitization and intracutaneous reactivity. All tests were performed in accordance with the International Standard Organization (ISO) 10993 Standard Series for Biological Evaluation of Medical Devices.

Invia – Silverlon NPWT Antimicrobial Wound Contact Dressing is not made with natural rubber latex.

Indications

Invia – Silverlon NPWT Antimicrobial Wound Contact Dressings are intended to provide moist wound environment and are used for the management of chronic, acute, surgical, traumatic, sub-acute and dehiscent wounds, ulcers (such as pressure or diabetic), partial-thickness burns, flaps and grafts. When used in conjunction with an NPWT pump that is intended to be used for negative pressure wound therapy (NPWT), the NPWT system is indicated for patients who benefit from wound management therapy via the application of negative pressure wound therapy for removal of fluids and excess exudates, infectious material, and tissue debris which may promote wound healing. It provides a non-adherent, conformable, antimicrobial barrier and may help reduce infection when placed under reticulated foam, gauze, or other wound fillers. Invia – Silverlon NPWT Antimicrobial Wound Contact Dressings may be cut or overlapped to fit any wound size.

Warnings and Precautions

Invia – Silverlon NPWT Antimicrobial Wound Contact Dressings are intended for external use only. Contact a health care professional if any of the following signs or symptoms are noted:

- increased pain, increased bleeding, increased swelling, increased wound drainage or increased redness in and around the wound site;
- there is a change in wound color and/or wound odor;
- the wound does not begin to show signs of healing; and
- any other unexpected symptoms occur.

- Consult a health care professional when Invia – Silverlon NPWT Antimicrobial Wound Contact Dressings are used with other wound care products.

- When using Invia – Silverlon NPWT Antimicrobial Wound Contact Dressings in conjunction with Invia negative pressure wound therapy systems, users should follow the directions for use supplied by the negative pressure wound therapy system manufacturer concerning the protection of exposed underlying structures.

- When using Invia – Silverlon NPWT Antimicrobial Wound Contact Dressing in conjunction with Invia negative pressure wound therapy systems, users should follow the directions for use supplied by the negative pressure wound therapy system manufacturer concerning specific wound conditions not suitable for negative pressure wound therapy and precautions to take in the event of the occurrence of bleeding during negative pressure wound therapy.

- Do not use past expiration date on the product packaging.

- Do not use if pouch is damaged or open.

- Do not use petroleum-based ointments or creams under Invia - Silverlon Dressings.

- Do not moisten Invia – Silverlon NPWT Antimicrobial Wound Contact Dressings with hydrogen peroxide or povidone iodine.

Some clinical studies have reported finding silver-resistant microbial strains when using silver based antimicrobial products. As of April, 2014 no adverse events or reports of Silverlon-resistant microbial strains have been received by Argentum.

Contraindications

- Do not use Invia – Silverlon NPWT Antimicrobial Wound Contact Dressings on patients with known sensitivity to silver or nylon.
- Do not use Invia – Silverlon NPWT - Antimicrobial Wound Contact Dressings on 3rd degree burns.

Adverse Reactions

N/A

Instructions for Use

- Cleanse wound with sterile water, distilled water, or normal saline, removing necrotic debris or eschar as needed per local protocol.
- Select the Invia – Silverlon NPWT Antimicrobial Wound Contact Dressing size that fully covers the wound up to the wound margins; trim the dressing edges as necessary.
- Activate Invia – Silverlon NPWT Antimicrobial Wound Contact Dressing by thoroughly moistening with sterile water, distilled water, or normal saline.
- Position Invia – Silverlon NPWT Antimicrobial Wound Contact Dressing directly over the wound, in contact with the skin; place Invia foam or gauze fill media of choice on top of the Invia – Silverlon NPWT Antimicrobial Wound Contact Dressing and secure in place according to the instructions supplied by the manufacture of the Invia negative pressure wound therapy equipment.
- Invia – Silverlon NPWT Antimicrobial Wound Contact Dressing may be used for three days, but may require frequent changing depending on wound condition and exudate build-up.
- Replace Invia – Silverlon NPWT Antimicrobial Wound Contact Dressing whenever the Invia foam / gauze fill media is replaced.
- To remove Invia – Silverlon NPWT Antimicrobial Wound Contact Dressing, first remove the foam / gauze fill media per local protocol, and then gently depress surrounding skin while lifting the edges of the Invia – Silverlon NPWT Antimicrobial Wound Contact Dressing.
 - If sticking of the dressing to the wound occurs, moisten the dressing as needed with sterile water, distilled water, or normal saline, until it can be easily removed by gently lifting the corners.

Storage

Store Invia – Silverlon NPWT Antimicrobial Wound Contact Dressings in normal warehouse conditions. Keep dry. Avoid excessive heat or humidity.

CAUTION: Federal U.S. law restricts this device to sale by or on the order of a physician.

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