



# Patient centered care without compromise

Leveraging our legacy of 60+ years of innovation in medical vacuum technology, our dynamic wound therapy solutions are versatile, dependable, and easy to use, ensuring continuity of care.



## Liberty and Motion Pumps



- Long battery life
- Compact size and lightweight design
- Whisper quiet

The Liberty and Motion NPWT pumps are billable under the CMS HCPCS code E2402.

Welcome to the world of  
Dynamic Wound Therapy

# Medaxis NPWT Dressings and Accessories

for use with Liberty and Motion NPWT pumps

Medaxis' entire portfolio of pumps and dressings is interchangeable with the Quick-connect, giving clinicians confidence they have a secure connection between pump and dressing.

## Integrated Dressing

Sterile "Peel and stick" for closed surgical incision and open wounds\*



Up to 7 days of wear time



Gentle silicone adhesive maintains secure seal



Evenly distributes pressure and draws off exudate to provide an optimal healing environment



## Foam and Gauze Dressing Kits



## Abdominal Dressing Kit

An effective solution<sup>1</sup> for the management of open abdominal wounds with exposed viscera and organs



## White Foam

Pre-moistened in sterile water with a higher tensile strength than black foam



## Y-Connector

Allows treatment of multiple wound sites on the same patient with only one NPWT device



## Silverlon®

Antimicrobial wound contact dressing



Caution: U.S. Federal law restricts this device to sale by or on the order of a physician



Medaxis AG  
Bahnhofstrasse 9  
6430 Baar  
Switzerland

Distributed by Medaxis LLC, 1101 Corporate Dr., McHenry, IL 60050 [medaxis.com](http://medaxis.com)  
1 Performance data on file at Medela AG  
\*Indicated for CSI and open wounds up to 0.5cm in depth without a wound filler  
Medela, Invia Motion, and Invia FitPad are registered with the US Patent and Trademark Office. Liberty is a trademark of Medela. Silverlon is a registered trademark of Argentum Medical, LLC and is used with permission.