

Wound healing that works for you

Medaxis delivers clinically effective, economically sustainable Negative Pressure Wound Therapy (NPWT) solutions that support health systems in achieving operational excellence, care continuity, and cost control across inpatient, outpatient, and home-based care settings.

Leveraging our legacy of 60+ years of innovation in medical vacuum technology, our dynamic wound therapy solutions are versatile, dependable, and easy to use, ensuring continuity of care.

Product Overview

Medaxis Liberty and **Motion NPWT** pumps incorporate our double-lumen Fitpad and electronically controlled feedback system to provide accurate pressure delivery and efficient fluid handling at the wound site.



Liberty NPWT Pump

Multi-user, reusable pump with a choice of pressure settings and therapy modes

-  2.2 lbs without canister
-  Air leak indicator
-  -40 to -200 mmHG
-  300 mL or 800 mL canister volume
-  Up to 14 hours
-  Audible alarms and warnings
-  Constant or intermittent



Motion NPWT Pump

Single-patient use pump that delivers traditional NPWT with a choice of pressure settings and therapy modes

-  0.88 lbs without canister
-  150 mL canister volume
-  -40 to -175 mmHG
-  Audible alarms and warnings
-  > 10 hours
-  Constant or intermittent

NPWT Dressings and Accessories

For use with Liberty and Motion NPWT pumps

Medaxis' entire portfolio of pumps and dressings is interchangeable with the Quick-connect, giving clinicians confidence they have a secure connection between pump and dressing.



Integrated Dressing

Sterile "peel-and-stick" dressing for closed surgical incisions and open wounds.

Up to 7 days of wear time.

Gentle silicone adhesive maintains a secure seal.

Evenly distributes pressure and draws off exudate to provide an optimal healing environment.



Foam and Gauze Dressing Kits

The foam and gauze are available in a kit with a FitPad and transparent film.

Foam and gauze can be cut to the appropriate size.



Abdominal Dressing Kit

Provides an effective solution¹ for the management of open abdominal wounds with exposed viscera and organs.



Y-Connector

Allows treatment of multiple wound sites on the same patient with only one NPWT pump.



White Foam

Pre-moistened in sterile water with a higher tensile strength than black foam.



Silverlon®

Antimicrobial wound contact dressing; provides a barrier to bacterial penetration using silver ions activated by moisture delivered in the dressing which can help reduce infection.



Caution: U.S. Federal law restricts this device to sale by or on the order of a physician

¹ Performance data on file at Medela AG. Medela, Invia Motion, and Invia FitPad are registered with the US Patent and Trademark Office. Silverlon is a registered trademark of Argentum Medical, LLC and is used with permission. 9000.3004_a Distributed by Medaxis LLC, 1101 Corporate Dr., McHenry, IL 60050